

BA

PatentWeb
HomeEdit
SearchReturn to
Patent List

Help

[Include in patent order

MicroPatent® Worldwide PatSearch: Record 1 of 1

[no drawing available]

Family Lookup**JP01118768****METHOD FOR STABILIZING COLOR FORMING TEST LIQUID**

WAKO PURE CHEM IND LTD

Inventor(s): ;SAKATA YOSHITSUGU ;HANADA TOSHIRO ;MUKAI TOYOJI**Application No.** 62276667 , **Filed** 19871031 , **Published** 19890511 ,

Abstract: PURPOSE: To enable more economical measurement with high accuracy by allowing cyclodextrin or/and the deriv. thereof to co-exist in a color forming test liquid contg. the oxidizable coloration reagent as a color forming component.

CONSTITUTION: The cyclodextrin or/and the deriv. thereof are allowed to co- exist in the color forming test liquid contg. the oxidizable coloration reagent expressed by the formula I or the formula II (where Z denotes an amino group, substd. amino group or hydroxy group; R¹ denotes a hydrogen atom, amino group, substd. amino group, hydroxy group, alkoxy group, alkyl group or alkenyl group; R²WR⁵ respectively independently denote a hydrogen atom, alkyl group, alkenyl group, acyl group, halogen atom, nitro group, sulfone group, carboxyl group, hydroxy group, or alkoxy group; R⁶ denotes a hydrogen atom, etc.; J denotes an oxygen atom, sulfur atom, etc.) as the color forming component. The trace components in a vital sample which can be measured by the color forming test liquid are exemplified by cholesterol, glucose, glycerol, free fatty acid, uric acid, bile acid, etc.

COPYRIGHT: (C)1989,JPO&Japio

Int'l Class: G01N03350 C12Q00126 C07D21906 C07D21908 C07D21916 C07D24146 C07D26538 C07D27930 C09B01500 C09B01700 C09B01900 C09B02100 D06P00148

MicroPatent Reference Number: 003070286

COPYRIGHT: (C) JPO

PatentWeb
HomeEdit
SearchReturn to
Patent List

Help

For further information, please contact:
Technical Support | Billing | Sales | General Information