Atty Dkt. No.: LIFE008

LFS-101

USSN: 09/593,827

a urea derivative dye on at least one surface of said matrix, wherein said urea derivative dye is 10-(carboxymethylaminocarbonyl)-3,7-bis(dimethylamino)phenothiazine or a salt thereof;

wherein said composition is stable for at least about six months at temperatures ranging from at least about -80°C to 60°C under humidity ranging from at least about 0% to 20%.

6. (Once Amended) The composition according to Claim 1, wherein said urea derivative dye is a member of a peroxide producing signal producing system present on said matrix.

11. (Twice Amended) A storage stable reagent test strip for use in detecting the presence or determining the concentration of an analyte in a physiological sample, said strip comprising:

a positively charged porous matrix comprising nylon; and

a peroxide producing signal producing system present on said matrix, wherein said peroxide producing signal producing system includes 10-(carboxymethylaminocarbonyl)-3,7-bis(dimethylamino)phenothiazine or a salt thereof,

ons (difficulty realisms) provided and a surface of the surface of

wherein said test strip is stable for at least about six months at temperatures ranging from at least about -80°C to 60°C under humidity ranging from at least about 0% to 20%.

- 19. (Twice Amended) An analyte detection or preasurement system comprising:
- (a) a storage stable reagent test strip comprising:
 - (i) a positively charged porous matrix comprising nylon; and
- (ii) a peroxide producing signal producing system present on said matrix, wherein said peroxide producing signal producing system includes 10-(carboxymethylaminocarbonyl)-3,7-bis(dimethylamino)phenothiazine or a salt thereof; and
- (b) an automated instrument,

wherein said test strip is stable for at least about six months at temperatures ranging from at least about -80°C to 60°C under humidity ranging from at least about 0% to 20%.

- 20. (Twice Amended) A method for detecting the presence or determining the concentration of an analyte in a sample, said method comprising:
 - (a) applying said physiological sample to a storage stable reagent test strip comprising:

Atty Dkt. No.: LIFE008

LFS-101

USSN: 09/593,827

(i) a positively charged porous matrix comprising aylon; and

(ii) a peroxide producing signal producing system present on said matrix, wherein said peroxide producing signal producing system includes 10-(carboxymethylaminocarbonyl)-3,7-bis(dimethylamino)phenothiazine or a salt thereof,

wherein said test strip is stable for at least about six months at temperatures ranging from at least about -80°C to 60°C under humidity ranging from at least about 0% to 20%;

- (b) detecting a signal produced by said signal producing system; and
- (c) relating said detected signal to the presence or concentration of said analyte in said physiological sample.
- 24. (Twice Amended) A kit for use in determining the concentration of an analyte in a physiological sample, said kit comprising:
 - (a) a storage stable reagent test strip comprising:
 - (i) a positively charged porous matrix comprising nylon; and
 - (ii) a peroxide producing signal producing system present on said matrix, wherein said peroxide producing signal producing system includes 10-(carboxymethylaminocarbonyl)-3,7-bis(dimethylamino)phenothiazine or a salt thereof,

wherein said test strip is stable for at least about six months at temperatures ranging from at least about -80°C to 60°C under humidity ranging from at least about 0% to 20%; and

- (b) at least one of;
 - (i) a means for obtaining said physiological sample and
 - (ii) an analyte standard.